

Flip Top Cap Centrifuge Tubes

50mL Centrifuge Tube, Flip Top Cap - Paperboard Rack, Sterile

Part Number: 229426

- ▶ Graduations are printed in bold, easy-to-read black ink and are accurate to +/-2%
- ▶ Large white marking area provides excellent contrast for identification
- ▶ Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- ▶ Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- ▶ Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- ▶ Ergonomic flip top cap is designed for single handed operation
- ▶ Paperboard rack is a recyclable, fold-flat alternative to bulky foam racks



Pack Description	25/Paperboard Rack with 1 Rack/Bag
Case Quantity	500
Sterile	Yes
Non-Pyrogenic	Yes
Autoclavable	Yes
Primary Color	Tube: Natural
Secondary Color	Cap: Green
Primary Material	Tube: Polypropylene
Secondary Material	Cap: HDPE
Maximum Volume	50mL
Graduated	Yes
Graduation Increments (Subdivisions)	2.5mL
Max Force (RCF)g	9400
DNase/RNase Free	Yes