

Erlenmeyer Flasks



Polycarbonate Erlenmeyer Flasks



PETG Erlenmeyer Flasks

- Designed for mammalian suspension culture applications as well as bacterial & yeast culture
- PETG available with either vent or solid cap to prevent common mistakes made with combination caps such as sealing the cap completely and preventing gas exchange
- Polycarbonate available with vent cap only
- Vent caps with over-sized vent area feature 0.22µm hydrophobic membrane which allows for improved gas exchange without contamination
- Avoid contamination with single-use, pre-sterilized flasks
- Lightweight material is easier to handle than glass
- Flasks include molded graduation lines
- Individually bagged for easy storage and handling
- Medical grade virgin Polycarbonate is autoclavable
- Medical grade virgin PETG reduces gas permeability

POLYCARBONATE ERLLENMEYER FLASKS - VENT CAP

A, S, R/D F, NP, PC

	Part No.	Product Description	Flask Material	Cap Type	Bottom Type	Max. Volume	Pack Description	Qty/ Case
NEW	229821	125mL Erlenmeyer Flask	Polycarbonate	Vent	Plain	125mL	1/Bag	12
NEW	229825	250mL Erlenmeyer Flask	Polycarbonate	Vent	Plain	265mL	1/Bag	12
NEW	229829	500mL Erlenmeyer Flask	Polycarbonate	Vent	Plain	525mL	1/Bag	12

PETG ERLLENMEYER FLASKS - VENT CAP

S, R/D F, NP, PETG

	229801	125mL Erlenmeyer Flask	PETG	Vent	Plain	125mL	1/Bag	12
	229805	250mL Erlenmeyer Flask	PETG	Vent	Plain	265mL	1/Bag	12
	229809	500mL Erlenmeyer Flask	PETG	Vent	Plain	525mL	1/Bag	12
	229813	1000mL Erlenmeyer Flask	PETG	Vent	Plain	1035mL	1/Bag	6

PETG ERLLENMEYER FLASKS - SOLID CAP

S, R/D F, NP, PETG

	229800	125mL Erlenmeyer Flask	PETG	Solid	Plain	125mL	1/Bag	12
	229804	250mL Erlenmeyer Flask	PETG	Solid	Plain	265mL	1/Bag	12
	229808	500mL Erlenmeyer Flask	PETG	Solid	Plain	525mL	1/Bag	12
	229812	1000mL Erlenmeyer Flask	PETG	Solid	Plain	1035mL	1/Bag	6