

# Snap-Pop Lid Centrifuge Tubes

## 50mL Centrifuge Tube, Snap-Pop Lid - Bag, Sterile

**Part Number:** 229498



- ▶ Single piece lid and tube available in conical and self-standing bottoms
- ▶ Tab on lid enables convenient one-handed opening
- ▶ Lid audibly snaps closed, and opens with a 'pop!' so you know when it is sealed
- ▶ Ideal for sample storage, collection, and experimental procedures

<b>Pack Description</b>	25/Re-sealable Bag
<b>Case Quantity</b>	500
<b>Sterile</b>	Yes
<b>Non-pyrogenic</b>	Yes
<b>Autoclavable</b>	Yes
<b>Primary Color</b>	Natural
<b>Primary Material</b>	Polypropylene
<b>Maximum Volume</b>	50mL
<b>Graduated</b>	Yes
<b>Graduation Increments (Subdivisions)</b>	2.5mL
<b>Max Force (RCF)g</b>	9000
<b>RNase/DNase Free</b>	Yes