

Snap-Pop Lid Centrifuge Tubes

50mL Centrifuge Tube, Snap-Pop Lid,
Self-Standing - Bag, Sterile

Part Number: 229499



- ▶ Single piece lid and tube available in conical and self-standing bottoms
- ▶ Tab on lid enables convenient one-handed opening
- ▶ Lid audibly snaps closed, and opens with a 'pop!' so you know when it is sealed
- ▶ Ideal for sample storage, collection, and experimental procedures
- ▶ Self-standing base for use without a rack

Pack Description	25/Re-sealable Bag
Case Quantity	500
Sterile	Yes
Non-pyrogenic	Yes
Autoclavable	Yes
Primary Color	Natural
Primary Material	Polypropylene
Maximum Volume	50mL
Graduated	Yes
Graduation Increments (Subdivisions)	2.5mL
Max Force (RCF)g	1000
RNase/DNase Free	Yes