

Flip Top Cap Centrifuge Tubes

15mL Centrifuge Tube, Flip Top Cap - Bag, Sterile

Part Number: 229417



- ▶ Graduations are printed in bold, easy-to-read black ink and are accurate to $\pm 2\%$
- ▶ Large white marking area provides excellent contrast for identification
- ▶ Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- ▶ Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- ▶ Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- ▶ Ergonomic flip top cap is designed for single handed operation

Pack Description	25/Re-sealable Bag
Case Quantity	500
Sterile	Yes
Non-Pyrogenic	Yes
Autoclavable	Yes
Primary Color	Tube: Natural
Secondary Color	Cap: Green
Primary Material	Tube: Polypropylene
Secondary Material	Cap: HDPE
Maximum Volume	15mL
Graduated	Yes
Graduation Increments (Subdivisions)	0.5mL
Max Force (RCF)g	9400
RNase/DNase Free	Yes