

# Centrifuge Tubes-Paperboard Rack

## 50mL Centrifuge Tube - Paperboard Rack, Sterile

Part Number: 229420



- ▶ Graduations are printed in bold, easy-to-read black ink and are accurate to +/-2%
- ▶ Large white marking area provides excellent contrast for identification
- ▶ Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- ▶ Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- ▶ Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- ▶ Paperboard rack is a recyclable, fold-flat alternative to bulky foam racks

<b>Pack Description</b>	25/Paperboard Rack with 1 Rack/Bag
<b>Case Quantity</b>	500
<b>Sterile</b>	Yes
<b>Non-Pyrogenic</b>	Yes
<b>Autoclavable</b>	Yes
<b>Primary Color</b>	Tube: Natural
<b>Secondary Color</b>	Cap: Green
<b>Primary Material</b>	Tube: Polypropylene
<b>Secondary Material</b>	Cap: HDPE
<b>Maximum Volume</b>	50mL
<b>Graduated</b>	Yes
<b>Graduation Increments (Subdivisions)</b>	2.5mL
<b>Max Force (RCF)g</b>	16000
<b>RNase/DNase Free</b>	Yes