

Centrifuge Tubes-Plastic Rack

15mL Centrifuge Tube - Plastic Rack, Sterile

Part Number: 229418

- ▶ Graduations are printed in bold, easy-to-read black ink and are accurate to +/-2%
- ▶ Large white marking area provides excellent contrast for identification
- ▶ Easy to write on flat surface polyethylene (HDPE) caps feature integrated leak-proof seal
- ▶ Tubes include graduations molded into the conical bottom, ideal for pellet viewing and small volume quantification
- ▶ Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- ▶ Plastic rack is manufactured from polypropylene and is autoclavable, reusable, and recyclable
- ▶ Chemically resistant to drips
- ▶ Design of bottom prevents racks from sticking/freezing to surfaces
- ▶ Widely spaced design catches drips and allows tubes to be retrieved more easily



Pack Description	25/Plastic Rack with 2 Racks/Bag
Case Quantity	300
Sterile	Yes
Non-Pyrogenic	Yes
Autoclavable	Yes
Primary Color	Tube: Natural
Secondary Color	Cap: Green
Primary Material	Tube: Polypropylene
Secondary Material	Cap: HDPE
Maximum Volume	15mL
Graduated	Yes
Graduation Increments (Subdivisions)	0.5mL
Max Force (RCF)g	12000
RNase/Dnase Free	Yes