

Rack Only

PLASTIC RACK ONLY, 15mL Centrifuge Tube, Reusable, Non-sterile

Part Number: 229419

- ▶ Plastic rack is manufactured from polypropylene and is autoclavable, reusable, and recyclable
- ▶ Chemically resistant to drips
- ▶ Design of bottom prevents racks from sticking/freezing to surfaces
- ▶ Widely spaced design catches drips and allows tubes to be retrieved more easily



Pack Description	5/Bag
Case Quantity	5
Sterile	No
Autoclavable	Yes
Primary Color	Green
Primary Material	Polypropylene