

# Rack Only

## PLASTIC RACK ONLY, 50mL Centrifuge Tube, Reusable, Non-sterile

**Part Number:** 229429

- ▶ Plastic rack is manufactured from polypropylene and is autoclavable, reusable, and recyclable
- ▶ Chemically resistant to drips
- ▶ Design of bottom prevents racks from sticking/freezing to surfaces
- ▶ Widely spaced design catches drips and allows tubes to be retrieved more easily



<b>Pack Description</b>	5/Bag
<b>Case Quantity</b>	5
<b>Sterile</b>	No
<b>Autoclavable</b>	Yes
<b>Primary Color</b>	Green
<b>Primary Material</b>	Polypropylene

# PLASTIC RACK ONLY, 50mL Centrifuge Tube, Reusable, 1

**Part Number:** 229429

- ▶ Plastic rack is manufactured from polypropylene and is autoclavable, reusable,
- ▶ Chemically resistant to drips
- ▶ Design of bottom prevents racks from sticking/freezing to surfaces
- ▶ Widely spaced design catches drips and allows tubes to be retrieved more easily



ted Products



**Autoclavable** 229421

**Primary Color**



**Serological Pipets**

Green



Part Number: 229429

Price: \$51.00

Qua

## Don't See What You're Looking For?

Let us help you find it! We frequently add new products and are committed to continuous improvement driven by customer feedback.

[CONTACT US](#)

### CELLTREAT Scientific Products

CELLTREAT Scientific Products is dedicated to manufacturing unique, high quality laboratory consumables at significant savings compared to alternative brands. Experience the C

#### Let Us Help You

[Quality Certificate](#)[Sample Request](#)[Distributors](#)[Resources](#)[Contact](#)[Employment](#)

#### Stay Connected

Sign up to stay informed of all the current products and new product updates.

Toll Free: 877.345.4620 | Phone: 978.650.1320 | Email: [info@ce](mailto:info@celltreat.com)



© 2020 Cytiva Scientific Products | All rights reserved

DOMINIQUE DUTSCHER SAS