

Spécification



SEFAR MEDIFAB 03-170/54 102 cm
Numéro d' article: 3A03-0170-102-10

Material	PA 6.6.
-----------------	---------

	Norme	Valeur	Tolerance	Unité
LAL		0	(+ 0,125 / - 0)	EU/ml
Hämo		0	(+ 0,03 / - 0)	hemolysis value
Laize		102	(+ 1 / - 1)	cm
Epaisseur	ISO 5084	100		µm
Ouverture de maille (chaîne)	ASTM E11	170	(+ 9 / - 9)	µm
Ouverture de maille (trame)	ASTM E11	170	(+ 9 / - 9)	µm
Nombre de fils (chaîne)	ASTM E11	43		1/cm
Nombre de fils (trame)	ASTM E11	43		1/cm
Diamètre du fil chaîne		62		µm
Diamètre du fil trame		62		µm
Surface libre		53,5		%
Poids		32		g/m2

Traitement	Medical washed
-------------------	----------------

Sefar provides industry leading, customized products mainly based on monofilament mesh and fabrics, expert services and professional customer support around the world for more than 180 years. Our products are used in a broad spectrum of filtration, separation and dosing applications.

The specifications of our products are designed to meet and exceed the challenging expectations of our customers.

Because of technological developments related to the products and the constant optimization of our manufacturing processes the specifications are subject to change without notice.

Please consult your local Sefar representative to verify that this information remains valid and check the suitability of this product for the targeted application before use.

Please note: This is an automatically generated document. All stated values without indications of tolerances are nominal values, which are not measured. All stated values with indications of tolerances are measured based on a product specific control plan and represent statistical mean values.