

Tissue Culture Test Plate 6 – 96 wells



Features of the Tissue Culture Test Plates:

- Lid with air-venting system for controlled gas exchange with low evaporation
- Vented base for air flow between stacked dishes to prevent condensation
- Sloped edge allows placement of the lid in one position only
- Yellow inscription field: yellow on yellow - match!
- Surface activated growth area on spherical zone of well only
- Crystal clear transparency for excellent viewing
- Black alphanumeric for quick identification of the wells
- Clear alphanumeric between wells for orientation during microscopy



Additional features with Tissue Culture Test Plates with 96 wells:

- Manufactured according to ANSI/SLAS 1-2004 to ANSI/SLAS 4-2004
- 2 well geometries available: F- and U-bottom.

The well geometry is an important criterion with a 96 well tissue culture plate. TPP offers 2 types of well geometries


F-base

- Excellent optical characteristics
- Suitable for precise optical measurements, the measuring light is not distracted by the geometry
- Suitable for microscopy applications (bottom reading)
- Growth area 0.322 cm²



U-base

- No rims, well suited for pipetting
- Used for agglutination tests and «+/-» evaluations
- Growth area: 0.636 cm²
- Working volume: 0.34 ml



1. Black alphanumerical labelling
2. Yellow inscription field
3. Excellent stacking
4. Alphanumerical identification mark next to each well
5. Gripping edges

Optimal gas exchange with lowest evaporation

Type	Product-No.	Number Wells pcs	Growth Area cm ²	Version	Volume ml	Inside- Ø mm	Dimensions l x w x h mm	Material	Qty/ Bag pcs	Qty/ Case pcs
	92006	6	9.026	F	15.53	33.9	128 x 86 x 22	PS	1	126
	92012	12	3.466	F	5.96	21.0	128 x 86 x 22	PS	1	126
	92024	24	1.864	F	3.18	15.4	128 x 86 x 22	PS	1	126
	92048	48	0.875	F	1.49	10.6	128 x 86 x 22	PS	1	126
	92096	96	0.322	F	0.36	6.4	128 x 86 x 17	PS	1	162
	92097	96	0.636	U	0.34	6.4	128 x 86 x 17	PS	1	162
	92406	6	9.026	F	15.53	33.9	128 x 86 x 22	PS	4	72
	92412	12	3.466	F	5.96	21.0	128 x 86 x 22	PS	4	72
	92424	24	1.864	F	3.18	15.4	128 x 86 x 22	PS	4	72
	92448	48	0.875	F	1.49	10.6	128 x 86 x 22	PS	4	72
	92696	96	0.322	F	0.36	6.4	128 x 86 x 17	PS	6	108
	92697	96	0.636	U	0.34	6.4	128 x 86 x 17	PS	6	108

Tissue Culture Test Plate 6 – 96 wells, single wrapped, small packaging unit

Type	Product-No.	Number Wells pcs	Growth Area cm ²	Version	Volume ml	Inside- Ø mm	Dimensions l x w x h mm	Material	Qty/ Bag pce	Qty/ Case pcs
	92106	6	9.026	F	15.53	33.9	128 x 86 x 22	PS	1	42
	92112	12	3.466	F	5.96	21.0	128 x 86 x 22	PS	1	42
	92124	24	1.864	F	3.18	15.4	128 x 86 x 22	PS	1	42
	92148	48	0.875	F	1.49	10.6	128 x 86 x 22	PS	1	42
	92196	96	0.322	F	0.36	6.4	128 x 86 x 17	PS	1	54
	92197	96	0.636	U	0.34	6.4	128 x 86 x 17	PS	1	54