


# CERTIFICATE OF CONFORMITY

**Product:** Shaker

**Tested by request of** N-BIOTEK, INC  
402-803 Technopark, 655 Pyeongcheon-ro, Wonmi-gu,  
Bucheon-si, Gyeonggi-do, Republic of Korea

**Manufactured at (name and place)** N-BIOTEK, INC  
402-803 Technopark, 655 Pyeongcheon-ro, Wonmi-gu,  
Bucheon-si, Gyeonggi-do, Republic of Korea

**Trade mark (if any):** 

**Model/Type Ref:** NB-101S RC

**Ratings and principal characteristics:** 230 V, 50 Hz, 17 W, Class I

**Additional information (if any) :**

A sample of the product has been tested and found to be in conformity with :

EN 61010-1:2010, EN 61010-2-051:2003

as shown in the test report(s) F690501/RF-SAF007211, F690501/RF-SAF007212

This Certificate of Conformity is the result of testing a sample of the product submitted, in accordance with the provisions of the relevant specific standard. It does not imply an assessment of the whole production. Conformity of the produced products with the specimen tested remains on the full responsibility of the manufacturer.

Brussels, 25/02/2015

  
Yuta Kim,  
Certification Manager.

