



Spex CertiPrep CLMS-1 Claritas PPT® Grade ICP-MS Multi-Element Solution 1, 10 µg/mL in 5% HNO₃; 125 mL

Spex CertiPrep – N° d'article FV-04568-80

Multi-element standards for ICP
and ICP-MS instrumentation

- Calibrate ICP-MS instrumentation prior to analysis
- Dilute to the concentration appropriate for the instrument
- Manufactured under ISO 9001, ISO/IEC 17025 and ISO 17034 compliance
- Multi-element internal standards in 5% HNO₃ matrix

Caractéristiques & Description

Volume (mL)	125
-------------	-----

COMPREND

Supplier Certificate of Analysis (CoA) is included with product shipment.