



Spex CertiPrep CLMS-4 Claritas PPT® Grade ICP-MS Multi-Element Solution 4, 10 µg/mL in H₂O/tr. HF/tr. HNO₃; 125 mL

Spex CertiPrep – N° d'article FV-04568-86

Multi-element standards for ICP
and ICP-MS instrumentation

- Calibrate ICP-MS instrumentation prior to analysis
- Designed to contain virtually every element in the mass spectrum for concentration verification checks
- Manufactured under ISO 9001, ISO/IEC 17025 and ISO 17034 compliance

Caractéristiques & Description

Volume (mL)	125
-------------	-----

COMPREND

Supplier Certificate of Analysis (CoA) is included with product shipment.