


TECHNICAL DATA SHEET

Part Number	999070
Description	Gel cutter tip



Benefits :

- ✓ Designed to remove bands and spots from agarose gels
- ✓ Avoids contamination
- ✓ One hand operation
- ✓ Ejector button tip release

Certificate :

- ✓ DNase and RNase free
- ✓ Can be autoclaved

Technical specifications :

- ✓ Fits standard 1 000 µl pipettes
- ✓ Dimensions : 1,1 x 6,5 mm

Sales unit :

250 Gel cutter tips