

Cell Culture water Pyrogen free

CAT N° : L0970

CAS number : 7732-18-5

Molecular weight : 18.02 g/mol

Chemical formula : H₂O

EC number : 231-791-2

Colour : colourless, clear solution

Storage conditions : Room temperature

Shelf life : 48 months

Sterility tests :

- Bacteria in aerobic and anaerobic conditions
- Fungi and yeasts

Endotoxin : < 0.005 EU/ml

Recommended use :

- Respect storage conditions of the product
- Do not use the product after its expiry date
- Store product in an area protected from light (not necessary for saline solutions).
- Manipulate the product in aseptic conditions (e.g. : under laminar air flow)
- Wear clothes adapted to the manipulation of the product to avoid contamination (e.g. : gloves, mask, hygiene cap, overall...)

The product is intended to be used in vitro, in laboratory only. Do not use it in therapy, human or veterinary applications.

Application :

This product is suitable to dissolve cell culture media in powder following the media technical data sheet.

Indications of deterioration :

This solution should be clear and free of particulate and flocculent material.

Do not use if the solution is cloudy or contains precipitate.

Other evidence of deterioration may include degradation of physical or performance characteristics.