

FICHE TECHNIQUE

Référence	NA934-1ML
Désignation	RABBIT IGG HRP LINKED WHOLE AB

Description :

- ✓ Amersham ECL Horseradish peroxidase (HRP) conjugated antibodies are highly species specific antibodies optimized for use with Amersham ECL Western blotting detection reagents.
- ✓ Highly species specific
- ✓ Optimized for use with our range of ECL Western blotting detection reagents
- ✓ Recommended dilutions minimize background
- ✓ Dilution guidelines are supplied for relevant applications. Store antibodies at 2 to 8°C.
- ✓ Application Western blotting
- ✓ Color Brown (light)
- ✓ Components Horseradish peroxidase conjugated antibody is supplied in phosphate buffered saline (sodium phosphate 0.1 M, NaCl 0.1 M) pH 7.5, containing 1% (w/v) bovine serum albumin and an anti-microbial agent
- ✓ Estimated shelf life from manufacture date 12 months
- ✓ The product is supplied with dilution guidelines for relevant applications
- ✓ Solubility : Soluble in cold water and hot water
- ✓ Source Donkey
- ✓ Storage conditions 2 - 8°C
- ✓ Wetted Part No
- ✓ Antibody form Whole Ab

Unité de vente :

1 ML rabbit igg hrp linked whole ab
