

Standard Analysis Certificate

PRODUCT :1,2-DICHLOROETHANE SPECTROSOL, for spectroscopy
CODE :P02827

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Refractive index at 20°C	-	1.443 - 1.447
Water content (K.F.)	mg/Kg	<= 100
Non volatile residue	mg/Kg	<= 5
Colour	Hazen	<= 10
Free acid (as HCl)	mg/Kg	<= 5
Assay (GC)	%	>= 99.8
UV transmittance at 225 nm	%	>= 10
UV transmittance at 230 nm	%	>= 50
UV transmittance at 240 nm	%	>= 90
UV transmittance at 250 nm	%	>= 95
UV transmittance at 260 nm	%	>= 99

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)