

## **$\alpha$ -tubulin Mouse Monoclonal Antibody**

### **OVERVIEW**

---

<b>CATALOG NO.</b>	PTX1032
<b>SIZE</b>	100 $\mu$ L
<b>HOST</b>	Mouse
<b>ISOTYPE</b>	IgG1
<b>CLONALITY</b>	Monoclonal Antibody

### **APPLICATIONS AND STORAGE**

---

<b>APPLICATIONS</b>	WB, IP, IF
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
<b>DELIVERY CONDITION</b>	Blue ice

**For research use only.**