

EBFP Mouse Polyclonal Antibody

OVERVIEW

CATALOG NO.	PTX1065
SIZE	100 µL
HOST	Mouse
ISOTYPE	IgG
CLONALITY	Polyclonal Antibody

APPLICATIONS AND STORAGE

APPLICATIONS	WB, IHC, IF, IP
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
DELIVERY CONDITION	Blue ice

For research use only.