

## **$\alpha$ tubulin Mouse Monoclonal Antibody(Zebrafish Specific)**

### **OVERVIEW**

---

|                    |                     |
|--------------------|---------------------|
| <b>CATALOG NO.</b> | PTX1078             |
| <b>SIZE</b>        | 100 $\mu$ L         |
| <b>HOST</b>        | Mouse               |
| <b>ISOTYPE</b>     | IgG1                |
| <b>CLONALITY</b>   | Monoclonal Antibody |

### **APPLICATIONS AND STORAGE**

---

|                                |  |
|--------------------------------|--|
| <b>APPLICATIONS</b>            | WB   |
| <b>STABILITY &amp; STORAGE</b> | PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.<br>4°C for short term (1 week), store at -20°C to -80°C for long term (1 year);<br>Avoid repeated freeze-thaw cycles |
| <b>DELIVERY CONDITION</b>      | Blue ice   |

**For research use only.**