

## GAPDH Mouse Monoclonal Antibody-HRP conjugated

### OVERVIEW

---

<b>CATALOG NO.</b>	PTX1088
<b>SIZE</b>	100 µL
<b>HOST</b>	Mouse
<b>ISOTYPE</b>	IgG1
<b>CLONALITY</b>	Monoclonal Antibody

### APPLICATIONS AND STORAGE

---

<b>APPLICATIONS</b>	WB, IHC
<b>STABILITY &amp; STORAGE</b>	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
<b>DELIVERY CONDITION</b>	Blue ice

For research use only.