

β -actin Mouse Monoclonal Antibody(Zebrafish Specific)

OVERVIEW

CATALOG NO.	PTX1095
SIZE	100 μ L
HOST	Mouse
ISOTYPE	IgG1
CLONALITY	Monoclonal Antibody

APPLICATIONS AND STORAGE

APPLICATIONS	WB, IHC, IF
STABILITY & STORAGE	PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. 4°C for short term (1 week), store at -20°C to -80°C for long term (1 year); Avoid repeated freeze-thaw cycles
DELIVERY CONDITION	Blue ice

For research use only.