

Human IgG1 Isotype Control antibody (5Vk2)

OVERVIEW

PRODUCT REFERENCE	PTX17956
CLONALITY	Monoclonal Antibody
ISOTYPE	IgG1

PRODUCT INFORMATION

EXPRESSION SYSTEM	Human
BUFFER	PBS pH7.4
STORAGE CONDITIONS	Store at 4 °C for short term. Freeze and store at -80 °C for long term. Addition of 20 to 40% glycerol may improve cryoprotection. Purified antibodies should ideally be stored as 0.5 to 2.0 mg/ml stock solutions. Aliquot the product according to your needs to avoid freeze/thaw cycles.