

# Standard Analysis Certificate

**PRODUCT** :ACETONE ATRASOL,for traces analysis  
**CODE** :P00532

TEST	U.M.	SPECIFICATION
Appearance	-	Clear colourless liquid
Refractive index at 20°C	-	1.357 - 1.361
Water content (K.F.)	mg/Kg	<= 500
Non volatile residue	mg/Kg	<= 2
Colour	Hazen	<= 10
Free acid (as CH <sub>3</sub> COOH)	mg/Kg	<= 20
Assay (GC)	%	>= 99.9
Ethanol	mg/Kg	<= 100
Methanol	mg/Kg	<= 500
2-Propanol	mg/Kg	<= 500
GC-ECD.Individual peak (CCl <sub>4</sub> )	µg/l	<= 1
Ret.range dichloromethane	-	-
to 1,2,4-trichlorobenzene	-	-
GC-ECD.Individual peak (Lindane)	ng/L	<= 2
Ret.range 1,2,4-trichlorobenzene	-	-
to decachlorobiphenyle	-	-
GC-FID.Individual peak (n-hexadecane)	µg/L	<= 2
Ret.range n-undecane to n-tetracontane	-	-

For traces analysis.

Date :11/01/2020

**QUALITY CONTROL RESPONSIBLE**  
**P. GUERARD (PEYPIN)**