

Standard Analysis Certificate

PRODUCT :DICHLOROMETHANE, stabilized with ethanol, ANHYDROUS ,analytical grade
CODE :P02910E

TEST	U.M.	SPECIFICATION
Refractive index at 20°C	-	1.422 - 1.426
Water content (K.F.)	mg/Kg	<= 50
Non volatile residue	mg/Kg	<= 10
Colour	Hazen	<= 10
Free acid (as HCl)	mg/Kg	<= 3
Assay (GC) (without stabilizer)	%	>= 99.95
Stabilizer (Ethanol)	% w/w	0.1 - 0.4
Carbon tetrachloride	mg/Kg	<= 20
Chloroform	mg/Kg	<= 40

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)