

Standard Analysis Certificate

PRODUCT :N,N-DIMETHYLFORMAMIDE For peptides synthesis
CODE :P03435

TEST	U.M.	SPECIFICATION
Appearance	-	Clear colourless liquid
Colour	Hazen	<= 10
Refractive index at 20°C	-	1.428 - 1.432
Water content (K.F.)	mg/Kg	<= 300
Bromophenol blue test	-	Conform
Amines content	mg/Kg	<= 5
Assay (GC)	%	>= 99.9
Methanol	mg/Kg	<= 100
Non volatile residue	mg/Kg	<= 15
Copper (Cu)	mg/kg	<= 0.05
Cadmium (Cd)	mg/Kg	<= 0.05
Chromium (Cr)	mg/Kg	<= 0.05
Iron (Fe)	mg/Kg	<= 0.05
Nickel (Ni)	mg/kg	<= 0.05
Lead (Pb)	mg/Kg	<= 0.1

DOMINIQUE DUTSCHER SAS

Date :11/01/2020

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)