

Notation produit / Product Rating

Référence / Part number	Désignation / Designation	Marque / Brand
P04-50500	HBSS 10x SANS CA, Mg, SANS NaHCO3 500ml	Pan Biotech

Type de produit / Product category

Consommable / Consumable

Fabrication / Production

Ce produit est fabriqué avec 30% d'énergie renouvelable. Pour les eaux WFI, Pan Biotech utilise le procédé d'osmose inverse, puis l'électrodéionisation et l'ultrafiltration au lieu de la distillation. Cela consomme moins d'énergie et s'avère plus efficace.

This product is manufactured with 30% renewable energy. The WFI production uses Reverse Osmosis followed by Electrodeionisation and Ultrafiltration instead of distillation. This process consumes less energy and is more efficient.

Emballage / Packaging

L'emballage de ce produit est en carton à poignées en fibres végétales. Pan Biotech propose de gros volumes de certains produits pour réduire le plastique (jusqu'à 70% de réduction pour 1L de produit) et récupère les bouteilles vides de milieux de culture pour recyclage.

The product packaging is made of cartons with handles made from grass fibres. Pan Biotech offers larger volumes of standard products to reduce the amount of PET (up to 70% for 1L). Empty PET cell culture media bottles can be returned to Pan Biotech for recycling.

Expédition / Shipping

Le produit est fabriqué en Allemagne.

The product is manufactured in Germany.

En savoir plus

[Le GoGreen chez Dutscher](#)

[Le développement durable chez Pan Biotech](#)