

## Datasheet

**L-Glutamine 200 mM****Cell Culture Supplement**

Product	Description	Catalogue-No.	Size
L-Glutamine	Cell culture supplement	P04-80050 P04-80100	50 ml 100 ml

**Product description**

L-Glutamine is an essential amino acid required by all mammalian and insect cells grown in culture. It serves as an additional source of carbon and energy and in particular is used for nucleic acid synthesis. Optimal cell performance usually requires supplementation of the medium with L-Glutamine prior to use.

L-Glutamine is more unstable in cell culture solutions than other amino acids. The rate and extent of L-glutamine degradation are related to storage temperatures, age of the product, and pH.

If you are unsure whether the concentration of L-Glutamine in your medium is enough add a small amount of L-Glutamine (~ 2 mM final concentration) to the culture medium. If the cell growth rate increases, L-Glutamine is most likely deficient and more should be added.

**Storage conditions**

Storage: -20° C

Stability: 2 years from date of production

**Composition**

L-Glutamine 200 mM in H<sub>2</sub>O

**Performance Testing**

L-Glutamine is tested for pH, osmolarity, and sterility. A Certificate of Analysis is available upon request.

**Technical support**

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email ([info@pan-biotech.com](mailto:info@pan-biotech.com)) or phone +49-8543-601630.

---

FOR RESEARCH USE ONLY! Not approved for human or animal diagnostic or therapeutic procedures.