

Standard Analysis Certificate

PRODUCT :TETRAHYDROFURAN ANHYDROUS ,analytical grade

CODE :P07010

TEST	U.M.	SPECIFICATION
Clear,colourless liq.appearance	-	Conform
Identification (IR)	-	Conform
Refractive index at 20°C	-	1.405 - 1.409
Density d20/4	-	0.884 - 0.894
Water content (K.F.)	mg/Kg	<= 100
Colour	Hazen	<= 10
Free acid (as CH ₃ COOH)	mg/Kg	<= 20
Peroxides (as H ₂ O ₂)	mg/Kg	<= 20
Stabilizer (ionol)	mg/Kg	200 - 400
Assay (GC)	%	>= 99.9
Non volatile residue (without stab.)	mg/Kg	<= 10

Stabilizer ionol : 2,6-di-tert-butyl-4-methylphenol ; BHT

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

P. GUERARD (PEYPIN)