

Standard Analysis Certificate

PRODUCT :TETRAHYDROFURAN ANHYDROUS ,over molecular sieve
CODE :P07010T

TEST	U.M.	SPECIFICATION
Clear, colourless liq. appearance	-	Conform
Identification (IR)	-	Conform
Refractive index at 20°C	-	1.405 - 1.409
Density d20/4	-	0.884 - 0.894
Water content (K.F.)	mg/Kg	<= 20
Colour	Hazen	<= 10
Free acid (as CH ₃ COOH)	mg/Kg	<= 10
Peroxides (as H ₂ O ₂)	mg/Kg	<= 20
Stabilizer (ionol)	mg/Kg	250 - 450
Assay (GC)	%	>= 99.7

Stabilizer ionol : 2,6-di-tert-butyl-4-methylphenol ; BHT - Over molecular sieve

Date

:11/01/2020

QUALITY CONTROL RESPONSIBLE

P. GUERARD (PEYPIN)