

Datasheet and Instructions for Use

**EDTA (1x)**

**EDTA in DPBS**

Product	Description	Catalogue-No.	Size
EDTA (1X)	EDTA 1:5000 in PBS	P10-15100	100 ml

**Product description**

EDTA (ethylenediaminetetraacetic acid) is a chelating agent. It is used to reduce cell clumping by removing Ca<sup>2+</sup> from cell wash and suspension media.

**Storage conditions**

Storage: 2 – 8 °C

Expiry: two years from date of production

**Composition**

0,2 g/L EDTA (Ethylenediaminetetraacetic acid)

PBS (w/o: Ca<sup>2+</sup> and Mg<sup>2+</sup>)

**Suitability**

FOR RESEARCH USE ONLY!

Not approved for human or animal diagnostic or therapeutic procedures.

**Technical Support**

Additional information will be available on our website: [www.pan-biotech.com](http://www.pan-biotech.com)

For technical support, questions or remarks please contact your local PAN-Biotech partner or the technical department of PAN-Biotech via email ([info@pan-biotech.com](mailto:info@pan-biotech.com)) or phone +49-8543-601630