

Standard Analysis Certificate

PRODUCT :SILICA 60 Å C.C 35-70 µm
CODE :P20000
METHOD :P20000

| TEST | U.M. | SPECIFICATION |
|-----------------------------|------|---------------|
| Particle size | µm | 35 ÷ 70 |
| Particle size <40 µm | % | <= 15 |
| Particle size >63 µm | % | <= 10 |
| Loss on drying | % | <= 8.0 |
| pH (5% <i>m/m</i> solution) | - | 6.5 ÷ 7.5 |

DOMINIQUE DUTSCHER SAS

Date :28/02/2015

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)