

Standard Analysis Certificate

PRODUCT :SILICA 60 A C.C 70-200 µm
CODE :P21000
METHOD :P21000

TEST	U.M.	SPECIFICATION
Particle size	µm	70 ÷ 200
Particle size <60 µm	%	<= 10
Particle size >200 µm	%	<= 10
Loss on drying	%	<= 6.0
pH (5%m/m solution)	-	6.5 - 7.5

DOMINIQUE DUTSCHER SAS

Date :28/02/2015

QUALITY CONTROL RESPONSIBLE
P. GUERARD (PEYPIN)