

# Standard Analysis Certificate

**PRODUCT** :1,3 - PROPANEDIOL PURE, Synthesis grade  
**CODE** :P80402

TEST	U.M.	SPECIFICATION
Appearance	-	Clear and viscous liquid
Identification (IR)	-	Conform
Refractive index at 20°C	-	1.438 - 1.442
Water content (K.F.)	mg/Kg	<= 1000
Colour	Hazen	<= 15
Assay (GC)	%	>= 99.7

DOMINIQUE DUTSCHER SAS

**Date** :11/01/2020

**QUALITY CONTROL RESPONSIBLE**  
**P. GUERARD (PEYPIN)**