

# Standard Analysis Certificate

**PRODUCT** :2-METHYL TETRAHYDROFURAN PURE, Synthesis grade  
**CODE** :P99602

TEST	U.M.	SPECIFICATION
Appearance	-	Clear colourless liquid
Refractive index at 20°C	-	1.404 - 1.408
Density d20/4	-	0.853 - 0.857
Water content (K.F.)	mg/Kg	<= 300
Peroxides (as H <sub>2</sub> O <sub>2</sub> )	mg/Kg	<= 100
Stabilizer (ionol)	mg/Kg	150 - 400
Assay (GC)	%	>= 99.9

Stabilizer ionol : 2,6-di-tert-butyl-4-methylphenol ; BHT

Date

:11/01/2020

**QUALITY CONTROL RESPONSIBLE**

**P. GUERARD (PEYPIN)**