

HyClone 100 mM Sodium Pyruvate solution



An easily accessible carbohydrate source for cell culture media.

An easily accessible carbohydrate source for cell culture media.

Part of a high-performance line of nutritional and growth supplements

Product Specifications [^]

Parameter	100 mM Sodium Pyruvate solution - 100 mL
Certifications/Compliance	Manufactured in cGMP (21 CFR 820) compliant and ISO9001 certified facilities
Endotoxin Concentration	≤ 10.00 EU/mL
Format	Liquid
Shelf life	12 months
Sterilization	0.1 µm sterile filtered
Pack size	100 mL
Order code	SH30239.01