

HyClone Bovine Serum Albumin



Cell culture media supplement.

Cell culture media supplement.

- pH 7.0
- Lyophilized powder

Product Specifications ^

SELECT YOUR MODEL

Bovine Serum Albumin (BSA) - 100 g

▼

Parameter	Bovine Serum Albumin (BSA) - 100 g
Certifications/Compliance	Packaged in cGMP (21 CFR 820) compliant and ISO9001 certified facilities
Endotoxin Concentration	< 10 EU/mg
Format	Powder
Sterilization	Non-sterile
Storage	2° to 8° C
pH	6.5 - 7.5
Order code	SH30574.02
Pack size	100 g