

## HyClone Amphotericin B (Fungizone) solution



An antifungal, anti-yeast and anti-mold reagent.

An antifungal, anti-yeast and anti-mold reagent.

- An antibiotic agent that is active against fungi
- Interferes with the permeability of the cell membrane of sensitive fungi by binding sterols
- Common working concentration is 2.5  $\mu\text{g}/\text{mL}$ , and it is cytotoxic at 30  $\mu\text{g}/\text{mL}$

Used as an antifungal, anti-yeast and anti-mold agent.

## Product Specifications ^

Amphotericin B (Fungiezone) - 50 mL



### Parameter

### Amphotericin B (Fungiezone) - 50 mL

Antibiotic Spectrum Antimycotic, Antifungal

Concentration 250 µg/mL

Format Liquid

Shelf life 18 months

Size 0.2 µm sterile filtered

Storage Store frozen

Pack size 50 mL

Order code SV30078.01