

## HyClone Antibiotic Antimycotic (Pen/Strep /Fungizone) solution



An antibiotic and antimycotic reagent.

An antibiotic and antimycotic reagent.

- Contains 10,000 units of penicillin, 10,000 µg of streptomycin, and 25 µg of Amphotericin B per mL
- Commonly used as an antifungal, anti-yeast, anti-mold or antibacterial
- Common working dose is 10 mL/L of 100X solution

Wide spectrum of activity against gram negative and gram positive bacteria as well as activity against fungi and yeast

## Product Specifications ^

---

Parameter	Antibiotic Antimycotic solution - 100 mL
Antibiotic Spectrum	Broad spectrum antibacterial and antifungal
Concentration	100X
Format	Liquid
Shelf life	18 months
Sterilization	0.2 µm sterile filtered
Storage	Store frozen
Pack size	100 mL
Order code	SV30079.01