

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tris-Borate-EDTA buffer - BioReagent, suitable for electrophoresis, 10× concentrate

Product Number: **T4415**

MDL: MFCD00236358



TEST	Specification
Appearance (Turbidity)	Clear
Appearance (Form)	Solution
pH	8.2 - 8.4
At 25 degrees Celcius	
Total Aerobic Microbial Count	≤ 10 CFU/ml
RNAse Detection	None Detected
When tested at a final concentration of 1X in the RNAT reaction	
Conc. by Titration (meq/g)	
Report result	
Suitability	Suitable
Suitable for use in Gel Electrophoresis after dilution to working concentration, and no DNase/Nickase activity detected.	
Expiration Date Period	Complete
12 months	

Specification: PRD.2.ZQ5.10000069480

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.