































## TECHNICAL DATASHEET

Part number	W14011
Description	Luer Valve Assortment Kit

		
Injection Site Male luer lock, Polycarbonate (PC), Qty:15	Check Valve, Polyester & Silicone, QTY: 5	Needleless Syringe-Activated Dual Check Valve, Polycarbonate & Silicone, Qty: 4
		
Needleless Syringe-Activated Check Valve, Polycarbonate, PVC & Silicone, Qty: 5	Slide Clamp for 1/8" Tubing, Polypropylene, Qty: 10	Roller Clamp for 3/16" OD Tubing, Polypropylene & Polystyrene, Qty: 5
		
Roller Clamp, Large Bore, Body PVC, Qty: 4	Pinch Clamp for 7mm Tubing, Polyacetal, Qty: 10	Female T Luer, Polycarbonate, Qty: 10
		
4-Port Infusion Y, Swivel Thread, Polycarbonate & PVC, Qty: 4	3-Port Infusion Y, Swivel Thread, Polycarbonate & ABS, Qty: 5	4-Way Stopcock, Polycarbonate & handle Polyethylene, Qty: 4
		
4-Way Luer Stopcock, Polycarbonate & Polyethylene, Qty: 4	4-Way Stopcock with Luer Lock, Polycarbonate & handle Polyethylene, Qty: 4	3-Way Stopcock, Polycarbonate & handle Polyethylene, Qty: 4
		
Pinch Clamp for 5mm Tubing, Polypropylene, Qty: 10	1-Way Stopcock, Polyacetal & Polycarbonate, Qty: 4	1-Way Stopcock, Polycarbonate & Polyethylene, Qty: 4
		
Manifold, 4-Port, Polycarbonate & Polyester Qty: 1	Manifold, 3-Port, Polycarbonate & Polyester, Qty: 1	Female Luer Fitting for 1/16" ID Tubing, Polypropylene, Qty: 20
		
Female Luer Fitting for 3/32" ID Tubing, Polypropylene, Qty: 20	Female Luer Fitting for 1/8" ID Tubing, Polypropylene, Qty: 20	Female Luer Fitting for 5/32" ID Tubing, Polypropylene, Qty: 20
		
Male Luer Fitting for 1/16" ID Tubing, Polypropylene Qty : 20	Male Luer Fitting for 3/32" ID Tubing, Polypropylene Qty: 20	Male Luer Fitting for 1/8" ID Tubing, Polypropylene Qty: 20
		
Male Luer Fitting for 1/8" ID Tubing, Polypropylene Qty: 20	Male/Female Luer Plug, Polycarbonate Qty: 20	Slide Clamp for 2.5mm O.D. Tubing, Polyethylene Qty: 20