Dutscher
ENSEMBLE, AIDONS LA SCIENCE À PROGRESSER

PRE-STACKED TAPED CASSETTES BIOSIGMA


Smart CHOICE
L'essentiel pour vos besoins, au meilleur prix.

Biosigma

Pre-stacked taped cassettes

Pre-stacked embedding cassettes arranged in strips of 40 pieces, held together by an easy-to-remove tape and designed for automated loading systems, ensuring consistent feeding, stable positioning, and reliable marking performance across common histology and biopsy procedures.

NEW



- Pre-stacked design with an easy-to-remove tape
- Lids included, supplied separately
- Stable positioning during automated loading
- 45-degree writing area for marking and printing
- Compatible with most stampers
- Resistant to histology solvents
- Available in multiple colors for efficient sample classification

Material :
- High-resistance POM (Polyoxymethylene)
- Chemical and thermal resistance: suitable for exposure to formalin, alcohol, xylene, and embedding temperatures



Packaging :
- Bottoms : 40 pcs/strip; 10 strips/box; 5 boxes/ctn
- Lids : 1000 pcs/box; 2 boxes/ctn

Use/Application :
- Suitable for routine processing, embedding, and preservation of tissues and biopsy samples in histology and pathology laboratories
- Suitable for use with automated or manual equipment, ensuring reliable handling throughout standard laboratory procedures

Model	Biopsy Cassettes with Tape	Histology Cassettes with Tape
White	390244	390249
Yellow	390245	390250
Pink	390246	390251
Green	390247	390252
Light blue	390248	390253
€ for 2000	NC -	NC -

Indicative prices, excluding VAT. Visit Dutscher.com

Dominique DUTSCHER SAS
2C, rue de Bruxelles - 67170 BERNOLSHEIM
Tél. 03 88 59 33 90 - Fax. 03 88 59 33 99
www.dutscher.com

Dutscher
ENSEMBLE, AIDONS LA SCIENCE À PROGRESSER