milian®
DUTSCHER GROUP

CENTRIFUGE TUBES

ClearLine®

Smart CHOICE


L'essentiel pour vos besoins, au meilleur prix.

From Routine to High-Speed Centrifugation !

- ClearLine® centrifuge tubes are designed to ensure reliable sample processing and consistent laboratory performance
- Manufactured from virgin transparent polypropylene and supplied with polyethylene screw caps featuring an internal sealing ring. Available as DailySpin tubes for routine centrifugation as well as High-Speed metal free tubes resistant up to 20,000 x g, for more demanding laboratory applications


ClearLine®
DailySpin Tubes

- For centrifugations up to 15,000 x g




ClearLine®
High Speed Metal-Free Tubes

- For centrifugations up to 20,000 x g



METAL FREE



Campus Spark
Chemin du Champ-des-Filles n°4A
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
milian@milian.com / www.milian.com

milian®
DUTSCHER GROUP


milian®
DUTSCHER GROUP

CENTRIFUGE TUBES

ClearLine®

Smart CHOICE
L'essentiel pour vos besoins, au meilleur prix.

ClearLine® DailySpin Tubes



Self-standing option with graduation at 2,5 ml

Volume
- 15 ml
- 50 ml

Material
- Tube : polypropylene
- Cap : polyethylene

Purity
- Human DNA / DNase / RNase / Endotoxin free / Non Cytotoxic

Sterilization
- Gamma sterilized (SAL 10^{-6})



Zip bag of 25 pieces



Rack of 25 pieces

Cat. No.	Description	Volume (ml)	Dimensions H x Ø (mm)	Packaging	Units/ box	CHF/ box
258372	Centrifuge conical tube in PP with screw cap	15	120 x Ø 16	Zip Bag of 25 pcs	500	NC -
258373	Centrifuge conical tube in PP with screw cap	15	120 x Ø 16	Rack of 25 pcs	500	NC -
258369	Centrifuge conical tube in PP with screw cap	50	116 x Ø 28	Zip Bag of 25 pcs	500	NC -
258370	Centrifuge conical tube in PP with screw cap	50	116 x Ø 28	Rack of 25 pcs	500	NC -
258371	Centrifuge self-standing tube in PP with screw cap	50	118 x Ø 28	Zip Bag of 25 pcs	500	NC -

Campus Spark
Chemin du Champ-des-Filles n°4A
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
milian@milian.com / www.milian.com

milian®
DUTSCHER GROUP

milian®
DUTSCHER GROUP


CENTRIFUGE TUBES

ClearLine®

Smart CHOICE
L'essentiel pour vos besoins, au meilleur prix.

ClearLine® High Speed Metal-Free Tubes

METAL FREE



Volume
- 15 ml
- 50 ml

Material
- Tube : polypropylene
- Cap : polyethylene


Purity
- Human DNA/DNase/RNase / Endotoxin free/Non Cytotoxic

Sterilization
- Gamma sterilized (SAL 10^{-6})

High-Speed Performance with Metal-Free Purity :
These centrifuge tubes are manufactured with a proprietary process to eliminate metal traces and prevent test interference.
Presence of elemental impurities is investigated through targeted analyses by Inductively-coupled Plasma/Mass Spectrometry (ICP-MS).

Metals and elements included in testing :
Aluminium, Antimony, Arsenic, Barium, Beryllium, Boron, Cadmium, Calcium, Chromium, Cobalt, Copper, Iron, Lead, Lithium, Magnesium, Manganese, Mercury, Molybdenum, Nickel, Palladium, Phosphorus, Platinum, Potassium, Rhodium, Selenium, Silver, Sodium, Strontium, Tellurium, Thallium, Tin, Titanium, Tungsten, Uranium, Vanadium, Zinc, Zirconium

METAL FREE


Zip Bag of 25 pieces


Rack of 25 pieces

Cat. No.	Description	Volume (ml)	Dimensions H x Ø (mm)	Packaging	Units/ box	CHF/ box
258455	Centrifuge conical tube in PP with screw cap	15	120 x Ø 28	Zip Bag of 25 pcs	500	NC -
258456	Centrifuge conical tube in PP with screw cap	15	120 x Ø 16	Rack of 25 pcs	500	NC -
258453	Centrifuge conical tube in PP with screw cap	50	116 x Ø 28	Zip Bag of 25 pcs	500	NC -
258454	Centrifuge conical tube in PP with screw cap	50	116 x Ø 28	Rack of 25 pcs	500	NC -

Find all ClearLine® products at milian.com

Campus Spark
Chemin du Champ-des-Filles n°4A
CH-1228 PLAN-LES-OUATES
Tél. +41 (0) 22 884 16 00
milian@milian.com / www.milian.com

milian®
DUTSCHER GROUP