Placas de Petri redondas

- De poliestireno
- Perfecta planitud
- Espesor constante
- Apilables
- Alta transparencia
- 2 versiones: no estéril (ISO 6) o estéril



Diámetro (mm)	55	60		90													150	
Altura (mm)	14	14,2		14,2		16,2		14,2		16,2		14,2		16,2		14,2	16,2	17,5
Lengüetas	3	3	Sin	3	3	3	3	Sin	Sin	Sin	Sin	1	1	1	1	3	3	3
Número de compartimentos								1								2	2	1
Estéril	Sí	Sí	Sí	No	Sí	No	Sí	No	Sí	No	Sí	No	Sí	No	Sí	No	Sí	Sí
Unid. de env.	10	10	10	20	20	20	20	20	20	20	20	20	20	20	20	20	20	5
Unid./caja	800	1080	1080	480	480	480	480	480	480	480	480	480	480	480	480	480	480	160
Ref.	390160*	390156	390152B	390155	390162	390153	390154	390163	390166	390164	390165	390167	390170	390168	390169	390158	390159	390157
€	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -

* Con rejilla de conteo.


Placa de Petri cuadrada

L x An (mm)	120 x 120
Altura (mm)	17,5
Lengüetas	3
Número de compartimentos	1
Estéril	Sí
Unid. de env.	10
Unid./caja	240
Ref.	390161
€	NC -



Placas de Petri de poliestireno


Diámetro (mm)	55		90		90		120		140	
Altura (mm)	14,2		14,2		16,2		15		15	
Lengüetas	Sin	3	Sin	3	Sin	3	Sin	3	Sin	3
Unid. de env.	15	15	25	25	25	25	16	16	10	16
Unid./caja	1005	1005	600	600	600	600	240	240	110	240
Ref. aséptica	483004	483005	483008	483009	483012	483013	483026	483027	483030	483031
€	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -
Ref. estéril	483006	483007	483010	483011	483014	483015	483028	483029	483032	483033
€	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -	NC -



Placas de Petri Rodac

- De poliestireno
- Con rejilla de conteo

Diámetro (mm)	65	
Altura (mm)	14,5	
Lengüetas	Sin	4
Número de compartimentos	1	
Unid. de env.	18	18
Unid./caja	918	918
Ref. aséptica	483000	483001
€	NC -	NC -
Ref. estéril	483002	483003
€	NC -	NC -



Placas de Petri específicas

Diámetro (mm)	90	90	90	90	90
Altura (mm)	14,2	14,2	14,2	14,2	14,2
Lengüetas	3	3	4	4	4
N° de compartimentos	1	2	1	2	3
Especificidad	Marcado láser		Anillo de estabilidad		
Unid. de env.	25	25	25	25	25
Unid./caja	600	600	600	600	600
Ref. aséptica	483016	483018	483020	483022	483024
€	NC -	NC -	NC -	NC -	NC -
Ref. estéril	483017	483019	483021	483023	483025
€	NC -	NC -	NC -	NC -	NC -



Precios indicativos. No incluyen IVA. Visiten nuestras WEB