
## Cubetas para espectrofotometría


NUEVO

### Window material transmission properties





| Materiales cubetas STARNA | Longitud de onda (nm) |
|---|---|
| OG = OPTICAL GLASS | 334 a 2500 nm |
| OS = SPECIAL OPTICAL GLASS | 320 a 2500 nm |
| QS = SPECTROSIL QUARTZ | 190 a 2700 nm |
| I = INFRASIL QUARTZ | 220 a 3800 nm |

Cubetas de espectrometría con 2 ventanas pulidas




| Paso óptico (mm) | Anl x L ext. (mm) | Al. con tapa /con tapón (mm) | An x L int. (mm) | Vol. (ml) | Ventanas pulidas | OG Ref. con tapa. PTFE | € | OG Ref. con tapón de PTFE | € | OS Ref. con tapa. PTFE | € | OS Ref. con tapón de PTFE | € |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| **Cubetas estándar** | | | | | | | | | | | | | |
| 1 | 12,5 x 3,5 | 45/55 | 10 x 1 | 0,4 | 2 | 244626 | NC - | 244578 | NC - | 244564 | NC - | 244559 | NC - |
| 2 | 12,5 x 4,5 | 45/55 | 10 x 2 | 0,7 | 2 | 244627 | NC - | 244579 | NC - | 244561 | NC - | 244586 | NC - |
| 5 | 12,5 x 7,5 | 45/48 | 10 x 5 | 1,7 | 2 | 244628 | NC - | 244580 | NC - | 244570 | NC - | 244534 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 10 x 10 | 3,5 | 2 | 244522 | NC - | 244581 | NC - | 244614 | NC - | 244531 | NC - |
| 20 | 12,5 x 22,5 | 45/48 | 10 x 20 | 7 | 2 | 244624 | NC - | 244582 | NC - | 244633 | NC - | 244587 | NC - |
| 40 | 12,5 x 42,5 | 45/48 | 10 x 40 | 14 | 2 | 244557 | NC - | 244583 | NC - | 244563 | NC - | 244588 | NC - |
| 50 | 12,5 x 52,5 | 45/48 | 10 x 50 | 17,5 | 2 | 244554 | NC - | 244584 | NC - | 244542 | NC - | 244555 | NC - |
| 100 | 12,5 x 102,5 | 45/48 | 10 x 100 | 35 | 2 | 244632 | NC - | 244585 | NC - | 244634 | NC - | 244589 | NC - |
| **Cubetas Semi-micro** | | | | | | | | | | | | | |
| 5 | 12,5 x 7,5 | 45/48 | 4 x 5 | 0,7 | 2 | - | - | - | - | - | - | 244637 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 4 x 10 | 1,4 | 2 | - | - | - | - | 244543 | NC - | 244560 | NC - |
| **Cubetas Micro** | | | | | | | | | | | | | |
| 5 | 12,5 x 7,5 | 45/48 | 2 x 5 | 0,35 | 2 | - | - | - | - | 244635 | NC - | - | - |
| 10 | 12,5 x 12,5 | 45/48 | 2 x 10 | 0,7 | 2 | - | - | - | - | 244636 | NC - | - | - |




| Paso óptico (mm) | Anl x L ext. (mm) | Al. con tapa /con tapón (mm) | An x L int. (mm) | Vol. (ml) | Ventanas pulidas | QS Ref. con tapa. PTFE | € | QS Ref. con tapón de PTFE | € | I Ref. con tapa. PTFE | € | I Ref. con tapón de PTFE | € |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| **Cubetas estándar** | | | | | | | | | | | | | |
| 1 | 12,5 x 3,5 | 45/55 | 10 x 1 | 0,4 | 2 | 244525 | NC - | 244524 | NC - | 244655 | NC - | 244594 | NC - |
| 2 | 12,5 x 4,5 | 45/55 | 10 x 2 | 0,7 | 2 | 244562 | NC - | 244533 | NC - | 244656 | NC - | 244595 | NC - |
| 5 | 12,5 x 7,5 | 45/48 | 10 x 5 | 1,7 | 2 | 244556 | NC - | 244540 | NC - | 244657 | NC - | 244596 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 10 x 10 | 3,5 | 2 | - | - | 244519 | NC - | 244658 | NC - | 244597 | NC - |
| 20 | 12,5 x 22,5 | 45/48 | 10 x 20 | 7 | 2 | 244625 | NC - | 244591 | NC - | 244659 | NC - | 244598 | NC - |
| 40 | 12,5 x 42,5 | 45/48 | 10 x 40 | 14 | 2 | 244638 | NC - | 244592 | NC - | 244660 | NC - | 244599 | NC - |
| 50 | 12,5 x 52,5 | 45/48 | 10 x 50 | 17,5 | 2 | 244553 | NC - | 244593 | NC - | 244661 | NC - | 244600 | NC - |
| 100 | 12,5 x 102,5 | 45/48 | 10 x 100 | 35 | 2 | 244567 | NC - | 244568 | NC - | 244662 | NC - | 244601 | NC - |
| **Cubetas Semi-micro** | | | | | | | | | | | | | |
| 10 | 12,5 x 12,5 | 45/48 | 4 x 10 | 1,4 | 2 | 244517 | NC - | 244532 | NC - | 244664 | NC - | 244700 | NC - |
| 40 | 12,5 x 42,5 | 45/48 | 4 x 40 | 5,6 | 2 | 244641 | NC - | 244683 | NC - | 244666 | NC - | 244702 | NC - |
| 100 | 12,5 x 102,5 | 45/48 | 4 x 100 | 14 | 2 | 244642 | NC - | | - | | - | | - |
| **Cubetas Semi-micro negras** | | | | | | | | | | | | | |
| 5 | 12,5 x 7,5 | 45/48 | 4 x 5 | 0,7 | 2 | 244643 | NC - | 244685 | NC - | - | - | - | - |
| 10 | 12,5 x 12,5 | 45/48 | 4 x 10 | 1,4 | 2 | 244527 | NC - | 244529 | NC - | 244668 | NC - | 244704 | NC - |
| 20 | 12,5 x 22,5 | 45/48 | 4 x 20 | 1,16 | 2 | 244644 | NC - | 244686 | NC - | - | - | - | - |
| 40 | 12,5 x 42,5 | 45/48 | 4 x 40 | 2,8 | 2 | 244645 | NC - | 244687 | NC - | - | - | - | - |
| 50 | 12,5 x 52,5 | 45/48 | 4 x 50 | 5,6 | 2 | 244646 | NC - | 244688 | NC - | - | - | - | - |
| **Cubetas Micro** | | | | | | | | | | | | | |
| 5 | 12,5 x 7,5 | 45/48 | 2 x 5 | 0,35 | 2 | 244647 | NC - | 244689 | NC - | 244669 | NC - | 244705 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 2 x 10 | 0,7 | 2 | 244526 | NC - | 244577 | NC - | 244670 | NC - | 244706 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 2 x 10 | 0,58 | 2 | 244530 | NC - | 244690 | NC - | 244671 | NC - | 244707 | NC - |
| 20 | 12,5 x 22,5 | 45/48 | 2 x 20 | 1,4 | 2 | 244648 | NC - | 244691 | NC - | 244672 | NC - | 244708 | NC - |
| 40 | 12,5 x 42,5 | 45/48 | 2 x 40 | 2,8 | 2 | 244649 | NC - | 244692 | NC - | 244673 | NC - | 244709 | NC - |
| 50 | 12,5 x 52,5 | 45/48 | 2 x 50 | 3,5 | 2 | 244650 | NC - | 244693 | NC - | 244674 | NC - | 244710 | NC - |
| **Cubetas Micro negras** | | | | | | | | | | | | | |
| 5 | 12,5 x 7,5 | 45/48 | 2 x 5 | 0,35 | 2 | 244651 | NC - | 244694 | NC - | 244675 | NC - | 244711 | NC - |
| 10 | 12,5 x 12,5 | 45/48 | 2 x 10 | 0,7 | 2 | 244576 | NC - | 244537 | NC - | 244676 | NC - | 244712 | NC - |
| 20 | 12,5 x 22,5 | 45/48 | 2 x 20 | 0,58 | 2 | 244548 | NC - | 244695 | NC - | 244677 | NC - | 244713 | NC - |
| 40 | 12,5 x 42,5 | 45/48 | 2 x 40 | 1,4 | 2 | 244652 | NC - | 244696 | NC - | 244678 | NC - | 244714 | NC - |
| 50 | 12,5 x 52,5 | 45/48 | 2 x 50 | 2,8 | 2 | 244653 | NC - | 244697 | NC - | 244679 | NC - | 244715 | NC - |
| 100 | 12,5 x 102,5 | 45/48 | 2 x 100 | 3,5 | 2 | 244654 | NC - | 244698 | NC - | 244680 | NC - | 244716 | NC - |

