
Cucharas espátulas







			De níquel		De acero inoxidable		De acero inoxidable, recubierto PTFE	
A (mm)	B (mm)	Vol. (ml)	Ref.	€	Ref.	€	Ref.	€
120	20	2	442195	NC -	442201	NC -	442205	NC -
150	20	2	442196	NC -	442202	NC -	442206	NC -
180	20	2	442197	NC -	442203	NC -	-	-
210	20	2	442198	NC -	442204	NC -	442207	NC -
100	8	0,2	442199	NC -	442208	NC -	442210	NC -
150	4	0,1	-	-	442209	NC -	-	-
150	8	0,2	442200	NC -	-	-	442211	NC -

Espátulas de acero inoxidable tipo análisis



Una de los extremos es rectangular y plana, y la otra en forma de cuchara. Dos longitudes: 235 mm y 500 mm.

Ref.	Longitud (mm)	€
001810	500	NC -

Espátulas de acero inoxidable




Espátulas diseñadas para el trasvase de polvo y materiales granulados.

Ref.	Descripción	Longitud (mm)	€
037872	Cuchara de acero inoxidable	225	NC -
037873	Espátula de acero inox. extremos redondo y plano	195	NC -
037874	Espátula, extremos redondo y plano	150	NC -

Espátulas tipo Chattaway





	Acero inoxidable		De acero inoxidable, recubierta PTFE	
L x An (mm)	Ref.	€	Ref.	€
100 x 8	442160	NC -	-	-
130 x 8	442161	NC -	-	-
150 x 8	442162	NC -	442163	NC -
200 x 8	442164	NC -	442165	NC -
100 x 10	442166	NC -	442174	NC -
130 x 10	442167	NC -	-	-
150 x 10	442168	NC -	442175	NC -
180 x 10	442169	NC -	442176	NC -
200 x 10	442170	NC -	442177	NC -
230 x 10	442171	NC -	442178	NC -
305 x 10	442172	NC -	-	-
350 x 10	442173	NC -	442179	NC -
100 x 4	442180	NC -	442185	NC -
130 x 4	442181	NC -	-	-
150 x 4	442182	NC -	442186	NC -
180 x 4	442183	NC -	-	-
200 x 4	442184	NC -	442187	NC -

Espátulas dobles flexibles



Espátula doble plana con dos extremos, flexible. Una espátula con borde redondeado y una con borde cuadrado. De acero inoxidable 18/8 de alta resistencia química y térmica, no magnética

Ref.	Longitud total (mm)	L x An espátula (mm)	€
076209	180	48 x 10	NC -
076210	210	48 x 10	NC -
076211	250	50 x 10	NC -
076212	210	52 x 20	NC -
076213	250	60 x 20	NC -
076214	300	65 x 22	NC -