Soportes inoxidables

Bases de acero inoxidable

- De acero inoxidable 18/10, DIN 12892
- Roscado M10
- Protección de caucho antideslizante

Ref.	L x An x Al (mm)	€
076573	210 x 130	NC -
076525	250 x 160	NC -
076574	300 x 150	NC -
076575	315 x 200	NC -

Bases de acero pintado

- De acero pintado DIN 12892
- Roscado M10

Ref.	L x An x Al (mm)	€
076576	180 x 100 x 6	NC -
076577	210 x 130 x 6	NC -
076578	250 x 160 x 6	NC -
076579	300 x 150 x 6	NC -
076580	315 x 200 x 6	NC -
076581	210 x 130 x 8	NC -
076582	250 x 160 x 8	NC -
076583	315 x 200 x 8	NC -
076584	315 x 200 x 10	NC -

Trípodes

- De fundición de acero pintado
- Roscado M10

Ref.	Longitud base (mm)	€
076585	105	NC -
076560	145	NC -

Nuez de acero inoxidable

- De acero inoxidable 18/10, DIN 12895
- Varilla de ø 16,5 mm



Ref.	Tipo	€
076603	1	NC -
076604	2	NC -
076605	3	NC -

Nuez de aluminio

- De aluminio, DIN 12895



Ref.	Tipo	€
076569	1	NC -

Pinzas especiales 3 dedos

- De acero inoxidable 18/10
- Dedos plastificados con silicona

Ref.	L x ø varilla (mm)	Apertura	€
076610	190 x 12	0 - 80	NC -
076611	190 x 12	0 - 120	NC -




076603


076569


076610



Varillas de acero inoxidable

- De acero inoxidable 18/10
- Para rosca M10

Ref.	ø x An (mm)	€
076521	12 x 500	NC -
076586	12 x 600	NC -
076568	12 x 750	NC -
076587	12 x 1000	NC -
076588	12 x 1250	NC -
076589	12 x 1500	NC -
076590	13 x 500	NC -
076591	13 x 600	NC -
076592	13 x 750	NC -
076593	13 x 1000	NC -
076594	13 x 1250	NC -
076595	13 x 1500	NC -
076596	16 x 1000	NC -
076597	16 x 1250	NC -
076598	16 x 1500	NC -

Varillas de acero zincado

- De acero zincado
- Para rosca M10

Ref.	ø x An (mm)	€
076599	12 x 500	NC -
076549	12 x 600	NC -
076561	12 x 750	NC -
076600	12 x 1000	NC -
076601	12 x 1250	NC -
076602	12 x 1500	NC -

Pinza de fijación



Ref.	Material	ø varilla (mm)	Tipo	€
076606	Acero inox.	16,5	1	NC -
076607	Acero inox.	16,5	3	NC -
076493	Aluminio	16,5	1	NC -
076608	Aluminio	16,5	3	NC -
076558	Aluminio	27	1	NC -

Precios indicativos. No incluyen IVA. Visiten nuestras WEB