Embudos Büchner de porcelana



Ref.	ø zona de filtración (mm)	ø agujeros (mm)	Altura total (mm)	€
140651	48	1	67	**NC -**
140652	62	1	94	**NC -**
140653	77	1	99	**NC -**
140654	97	1,5	111	**NC -**
140655	116	1,5	132	**NC -**
140656	130	2	137	**NC -**
140657	156	2	155	**NC -**

Embudos Büchner de polipropileno



- De polipropileno autoclavable
- Desmontables en dos partes para una limpieza más fácil



Diámetro int. sup. (mm)	59	76	93	130	172
Diámetro ext. varilla (mm)	10	11	12	22	20
Longitud varilla (mm)	45	65	71	93	107
Altura total (mm)	104	146	169	213	283
Capacidad (ml)	83	186	358	919	2906
Diámetro papel filtro (mm)	55	70	90	110	150
Unid./ bolsa	1	1	1	1	1
Ref.	**029511**	**029512**	**029513**	**029514**	**029515**
€/unidad	**NC -**	**NC -**	**NC -**	**NC -**	**NC -**

Embudos Büchner de polipropileno



- Autoclavables
- Moldeados en 2 piezas que se pueden separar fácilmente
- Resistentes a la corrosión y al calor
- Excelente resistencia química



Ref.	Volumen (ml)	Diámetro (mm)	Longitud (mm)	Diámetro de los agujeros (mm)	€
391586	40	42,5	47	1,2	**NC -**
391587	70	55	57	1,1	**NC -**
391471	180	70	57	2	**NC -**
391588	285	80	65	2	**NC -**
391589	390	90	68	2,5	**NC -**
391590	810	110	92	2,5	**NC -**
391591	2100	160	105	2,75	**NC -**
391592	6000	240	143	3	**NC -**

Crisoles filtrantes de vidrio Duran

DWK LIFE SCIENCES



- Vidrio de borosilicato 3,3
- Para el secado de precipitados
- Para completar con junta cónica 046571

Porosidad			100 - 160 µm		40 - 100 µm		16 - 40 µm		10 - 16 µm		1 - 1,6 µm	
Volumen (ml)	Ø (mm)	Unid./ caja	Ref.	€/ caja	Ref.	€/ caja	Ref.	€/ caja	Ref.	€/ caja	Ref.	€/ caja
8	24	10	-	-	**083450**	**NC -**	**084346**	**NC -**	**084347**	**NC -**	-	-
15	28	10	**084348**	**NC -**	**084349**	**NC -**	**084350**	**NC -**	**084351**	**NC -**	-	-
30	36	10	**084352**	**NC -**	**084353**	**NC -**	**092103**	**NC -**	**084354**	**NC -**	**084355**	**NC -**
50	46	10	**084356**	**NC -**	**084357**	**NC -**	**092104**	**NC -**	**092117**	**NC -**	**084358**	**NC -**



Juego de juntas de goma cónicas para embudos de filtración. Piezas de 22 a 84 mm.

Ref.	Descripción	€
046571	Juego de 7 juntas	**NC -**