Matraces aforados Hirschmann

A - Matraces aforados clase A, tapón poly graduación marrón



Ref.	Volumen (ml)	Tipo	Tolerancia (+/-ml)	Altura (+/- 5 mm)	Unid./ caja	€/ caja
501031	5	NS 7/16	0,025	70	2	NC -
501032	5	NS 10/19*	0,04	70	2	NC -
501033	10	NS 7/16	0,025	90	2	NC -
501034	10	NS 10/19*	0,04	90	2	NC -
501035	20	NS 10/19	0,04	110	2	NC -
501036	20	NS 12/21	0,06	105	2	NC -
501037	25	NS 10/19	0,04	110	2	NC -
501038	25	NS 12/21*	0,06	110	2	NC -
501039	50	NS 12/21	0,06	140	2	NC -
501040	50	NS 14/23*	0,08	140	2	NC -
501041	100	NS 12/21	0,1	170	2	NC -
501042	100	NS 14/23	0,1	170	2	NC -
501043	200	NS 14/23	0,15	210	2	NC -
501044	250	NS 14/23	0,15	220	2	NC -
501045	500	NS 19/26	0,25	260	2	NC -
501046	1000	NS 29/32*	0,6	300	1	NC -
501047	1000	NS24/29	0,4	300	1	NC -
501048	2000	NS29/32	0,6	370	1	NC -
501049	5000	NS34/35	1,2	470	1	NC -
501050	10000	NS45/40	2,0	570	1	NC -

* Boca ancha.

E - Matraces aforados vidrio marrón clase A

- Certificado de conformidad (con número de lote fechado)
- Con tapón de polipropileno o vidrio
- Conforme normas DIN EN ISO 1042



					Frasco con tapón de plástico		Frasco tapón de vidrio	
Capacidad (ml)	Tipo	Tolerancia (+/-ml)	Altura (+/- 5 mm)	Unid./ caja	Ref.	€/ caja	Ref.	€/ caja
5	NS 7/16	0,025	70	2	500988	NC -	501003	NC -
5	NS 10/19*	0,04	70	2	500989	NC -	501004	NC -
10	NS 7/16	0,025	90	2	500990	NC -	501005	NC -
10	NS 10/19*	0,04	90	2	500991	NC -	501006	NC -
20	NS 10/19	0,04	110	2	500992	NC -	501007	NC -
25	NS 10/19	0,04	110	2	500993	NC -	501008	NC -
50	NS 12/21	0,06	140	2	500994	NC -	501009	NC -
50	NS 14/23*	0,08	140	2	500995	NC -	501010	NC -
100	NS 12/21	0,1	170	2	500996	NC -	501011	NC -
100	NS 14/23	0,1	170	2	500997	NC -	501012	NC -
200	NS 14/23	0,15	210	2	500998	NC -	501013	NC -
250	NS 14/23	0,15	220	2	500999	NC -	501014	NC -
500	NS 19/26	0,25	260	2	501000	NC -	501015	NC -
1000	NS 24/29	0,4	300	1	501001	NC -	501016	NC -
2000	NS 29/32	0,6	370	1	501002	NC -	501017	NC -

* Boca ancha.

F - Matraces aforados clase A, vidrio marrón USP

- De vidrio borosilicato marrón
- Con tapón de polipropileno



Ref.	Volumen (ml)	Tolerancia (+/-ml)	Altura (+/- 5 mm)	Unid./ caja	€/ caja
501071	5	0,02	70	2	NC -
501072	10	0,02	90	2	NC -
501073	20	0,03	110	2	NC -
501074	25	0,03	110	2	NC -
501075	50	0,05	140	2	NC -
501076	100	0,08	170	2	NC -
501077	200	0,1	210	2	NC -
501078	250	0,12	220	2	NC -
501079	500	0,15	260	2	NC -
501080	1000	0,3	300	1	NC -
501081	2000	0,5	370	1	NC -

D - Matraces aforados de clase A, tapón trapezoidal

- Certificado de conformidad (con número de lote fechado)
- Con tapón de polipropileno o vidrio forma trapecio
- Conforme normas DIN EN ISO 1042



				Tapón polipropileno		Tapón vidrio	
Vol. (ml)	Tolerancia (+/-ml)	Tolerancia (+/-ml)	Altura (+/- 5 mm)	Ref.	€ 2 unid.	Ref.	€ 2 unid.
1	NS 7/16	0,025	60	501051	NC -	501060	NC -
2	NS 7/16	0,025	70	501052	NC -	501061	NC -
5	NS 7/16	0,025	70	501053	NC -	501062	NC -
5	NS 10/19*	0,04	70	501054	NC -	501063	NC -
10	NS 7/16	0,025	90	501055	NC -	501064	NC -
10	NS 10/19*	0,04	90	501056	NC -	501065	NC -
20	NS 10/19	0,04	110	501057	NC -	501066	NC -
25	NS 10/19	0,04	110	501058	NC -	501067	NC -
50	NS 12/21	0,06	140	501059	NC -	501068	NC -

* Boca ancha.

Precios indicativos. No incluyen IVA. Visiten nuestras WEB